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The screenshot displays the NuGenesis ELN software interface. The main window shows a document titled "試験記録書 (溶出)" (Test Record Book (Elution)). The left sidebar contains a table of contents with items 1 through 9, where "6. HPLC結果" (HPLC Results) is selected. The main area displays the "HPLC結果" form, which includes a header section with fields for "プロジェクト名" (Project Name), "結果セットビュー表示" (Result Set View Display), "結果セットID" (Result Set ID), "サンプルセット開始日" (Sample Set Start Date), and "サンプルセット終了日" (Sample Set End Date). Below this is a table for "結果取得" (Result Acquisition) with columns for "サンプル名" (Sample Name), "結果ID" (Result ID), "面積 (μV*秒)" (Area (μV*sec)), "バイアル" (Vial), "注入量 (μL)" (Injection Volume (μL)), "分析時間 (分)" (Analysis Time (min)), "分析日" (Analysis Date), "装置メソッドID" (Instrument Method ID), and "解析メソッドID" (Analysis Method ID). The table has 14 rows. The right sidebar shows the "ドキュメント詳細" (Document Details) section, which includes fields for "ドキュメントID" (Document ID), "ステータス" (Status), "タイトル" (Title), "サブタイトル" (Subtitle), "説明" (Description), "作成日" (Creation Date), "最終変更" (Last Modified), "デパートメント" (Department), "チーム" (Team), "必要なトレーニング" (Required Training), "テンプレート" (Template), "テンプレートカテゴリ" (Template Category), "設定" (Settings), "フォーム" (Form), "連番" (Serial Number), "印刷コントロール" (Print Control), "印刷履歴" (Print History), "デフォルト署名者" (Default Signer), "デフォルト確認者" (Default Checker), and "デフォルト承認者" (Default Approver).

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The screenshot shows the NuGenesis ELN software interface. The main window displays a table of sample information. A yellow callout box with the text "現物バーコードラベルを読み取り、サンプルデータを登録。" (Read the actual barcode label and register sample data.) points to a barcode reader icon. The table shows two samples with their respective lot numbers, item numbers, types, weights, units, and sample IDs.

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Projects

- Links
- Recently viewed
- All Available
 - NUGENESIS
 - AQT_Reports
 - AQT_Sandbox
 - Archived Empower Projects
 - Empower_Reports
 - Mass Spec Files
 - Mass Spec Reports
 - PrintCapture_Demo
 - SDMS_Administration
 - SDMS_System_Activities
 - Systems QT Projects**
 - All Projects
 - Archived Last Month
 - Archived This Month
 - Archived Today
 - DefaultView
 - Past 180 Days
 - Past 30 Days

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
Project: AQT_Reports View: All Reports


Select All None FILTER APPLIED 0

	ID	Empower Node	User Name	System Name	Report Name	Project Name	Emp DB Service Name	Sample
<input checked="" type="checkbox"/>	1	237	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna
<input type="checkbox"/>	2	232	Wli2k3vm	Demonstration User (demouser)	Alliance	LC Calibration Report	WLI2K3VM.WORLD	Externa
<input type="checkbox"/>	3	230	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna
<input type="checkbox"/>	4	225	Wli2k3vm	Demonstration User (demouser)	Alliance	LC Calibration Report	WLI2K3VM.WORLD	Externa
<input type="checkbox"/>	5	220	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna
<input type="checkbox"/>	6	215	Wli2k3vm	Demonstration User (demouser)	Alliance	LC Calibration Report	WLI2K3VM.WORLD	Externa
<input type="checkbox"/>	7	213	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna
<input type="checkbox"/>	8	208	Wli2k3vm	Demonstration User (demouser)	Alliance	LC Calibration Report	WLI2K3VM.WORLD	Externa
<input type="checkbox"/>	9	202	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna
<input type="checkbox"/>	10	197	Wli2k3vm	Demonstration User (demouser)	Alliance	LC Calibration Report	WLI2K3VM.WORLD	Externa
<input type="checkbox"/>	11	195	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna
<input type="checkbox"/>	12	190	Wli2k3vm	Demonstration User (demouser)	Alliance	LC Calibration Report	WLI2K3VM.WORLD	Externa
<input type="checkbox"/>	187	Wli2k3vm	System	486_sys	LC Calibration Report	WLI2K3VM.WORLD	Interna	

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NuGenesis VISION Preview - [NUGENESIS -- AQT_Reports]

File Edit View Send Options Annotations Launch Help

ID: 237
 User Name: Demonstration User (demonstrator)
 Report Name: LC Calibration Report
 Emp DB Service Name: WLI2K3VM.WORLD
 Integration Algorithm: ApexTrack
 Date Printed (Report): Wednesday, May 04, 2011
 Sample Name: NA
 Study Description:
 Processing Method: Internal Std Processing

Empower Node: Wli2k3vm
 System Name: 486_sys
 Project Name:
 Sample Set Name: Internal Std Sample Set
 Date/Time of Print: 2011-05-04 14:48:10
 Date Acquired: NA
 Study Name:
 Acquired Method Set: Internal Std Method Set
 Build Version: Empower 2 Software Build 2154

Date Calibrated: Wednesday, May 04, 2011 2:48:09 PM EDT

Amount

Area

Name: vanillin; RT: 6.500; Fit Type: Linear (1st Order); Cal Curve Id: 2335; R: 0.999099; R^2: 0.999817; Weighting: None; Equation: Y = 9.18e-001 X - 2.00e-001; Calibration Id 2334

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
All Reports in AQT_Reports on NUGENESIS - NuGenesis VISION - Windows Internet Explorer

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☐ Find any keyword (OR)
☒ Find all keywords (AND)

Current Selected All

Project: AQT_Reports View: All Reports

Select All None SEARCH APPLIED 0

	ID	Empower Node	User Name	System Name	Report Name
<input type="checkbox"/> 1	238	Wli2k3vm	Demonstration User (demouser)	486_sys	Summary Report for Instd OQ
<input type="checkbox"/> 2	237	Wli2k3vm	Demonstration User (demouser)	486_sys	LC Calibration Report
<input type="checkbox"/> 3			Demonstration User (demouser)	486_sys	Summary Report for Instd OQ
<input type="checkbox"/> 4			Demonstration User (demouser)	486_sys	LC Calibration Report
<input type="checkbox"/> 5			Demonstration User (demouser)	486_sys	Summary Report for Instd OQ
<input type="checkbox"/> 6			Demonstration User (demouser)	486_sys	LC Calibration Report
<input type="checkbox"/> 7	188	Wli2k3vm		486_sys	Summary Report for Instd OQ
<input type="checkbox"/> 8	187				
<input type="checkbox"/> 9	173	Wli2k3vm	System	486_sys	Summary Report for Instd OQ

All Content
Record Content
Record Tags
eSig Data
Annotation Data
Extended Details

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Select the Data View of the audit information

- Archive Device
- Archive Device
- Archive Exception Lists
- Archive Templates
- Results**
- Email Notifications
- Legal Holds
- Notifier List
- Policies
- Privileges
- Projects
- Reason
- Script
- Search Groups
- Servers
- Templates
- User Groups
- Users
- Views
- Volume Life Cycle

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Versions

Version Number	Version Type	Reason	Modified on	Modification
1	M		09.05.2013 08:38:22 by Tadahiro Kurita	Created From template
1	M		09.05.2013 09:05:03 by Tadahiro Kurita	Updated
1	R		09.05.2013 09:05:11 by Tadahiro Kurita	Released
1	A	OK	09.05.2013 09:05:24 by Tadahiro Kurita	Approved
1	H	OK	09.05.2013 09:05:51 by Tadahiro Kurita	History
2	M		09.05.2013 09:05:51 by Tadahiro Kurita	Created
2	R		09.05.2013 09:05:56 by Tadahiro Kurita	Released

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
Hide Method configuration Highlight input fields Switch to design/run-time mode Close Laboratory Execution Method Editor

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IPM Prop Inventory

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Laboratory Execution Method Editor

- Format Form
Format the section
- Sub Header
add a subheader
- Input field
add an input
- Procedure
add a procedure
- Form Settings
General Form Settings
- Cell Settings
Cell specific settings as metadata transfer, validation, lists, SDS operations
- Library
Import/Export an entire form configuration

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Copyright (C) Business Engineering Corporation.